

APPENDIX III

GEL T1 AND T2 REFERENCE VALUES

for various

FIELD STRENGTHS

and

TEMPERATURES

EUROSPIN II MR TEST OBJECTS

AGAROSE GEL SAMPLES

The Gel Samples enclosed in the Tubes with TO5 are made from Agarose doped with Gadolinium. This represents an ideal test material for MR Imaging Studies since it is possible, by varying the composition, to produce gels of any desired T1 and T2.

Additionally, their method of preparation is so well controlled that a mathematical model has been developed of the relaxation process. Thus it is possible to predict the T1 and T2 at all frequencies and temperatures within the imaging range.

The data presented in this section comes from this model, but each gel has been subject to rigid quality control and the relaxation times have been checked at a spot frequency (8MHz).

If values are required for fields other than those stated (and temperatures) please contact Diagnostic Sonar Ltd for information.

APPENDIX III
GEL T1 AND T2 VALUES

Imaging field 0.02 Tesla

| Temperature (K) | | | | | | |
|-----------------|---------|-----|------|---------|-----|-----|
| | 292 | 296 | 300 | 292 | 296 | 300 |
| Tube No | T1 (ms) | | | T2 (ms) | | |
| 1 | 87 | 92 | 97 | 41 | 40 | 39 |
| 2 | 132 | 140 | 148 | 59 | 56 | 55 |
| 3 | 137 | 145 | 154 | 83 | 82 | 81 |
| 4 | 185 | 193 | 202 | 47 | 44 | 42 |
| 5 | 201 | 213 | 224 | 79 | 75 | 72 |
| 6 | 211 | 224 | 238 | 117 | 116 | 114 |
| 7 | 267 | 281 | 295 | 83 | 78 | 75 |
| 8 | 280 | 296 | 313 | 114 | 109 | 105 |
| 9 | 343 | 361 | 380 | 103 | 97 | 92 |
| 10 | 358 | 379 | 401 | 134 | 128 | 123 |
| 11 | 425 | 448 | 471 | 123 | 115 | 110 |
| 12 | 883 | 944 | 1008 | 347 | 333 | 322 |
| 13 | 500 | 530 | 562 | 192 | 184 | 177 |
| 14 | 510 | 538 | 566 | 150 | 141 | 134 |
| 15 | 607 | 643 | 680 | 192 | 181 | 173 |
| 16 | 668 | 706 | 744 | 186 | 174 | 166 |
| 17 | 644 | 679 | 714 | 168 | 157 | 149 |
| 18 | 719 | 756 | 791 | 162 | 151 | 142 |

Estimated accuracy = $\pm 3\%$

APPENDIX III

GEL T1 AND T2 VALUES

Imaging field 0.15 Tesla

| Temperature (K) | | | | | | |
|-----------------|---------|------|------|---------|-----|-----|
| | 292 | 296 | 300 | 292 | 296 | 300 |
| Tube No | T1 (ms) | | | T2 (ms) | | |
| 1 | 118 | 124 | 130 | 45 | 43 | 42 |
| 2 | 180 | 189 | 198 | 65 | 61 | 59 |
| 3 | 180 | 189 | 198 | 94 | 92 | 90 |
| 4 | 275 | 289 | 304 | 50 | 46 | 43 |
| 5 | 276 | 290 | 306 | 86 | 81 | 77 |
| 6 | 277 | 291 | 307 | 132 | 128 | 125 |
| 7 | 376 | 397 | 419 | 89 | 83 | 78 |
| 8 | 378 | 399 | 421 | 123 | 117 | 112 |
| 9 | 483 | 510 | 540 | 109 | 102 | 96 |
| 10 | 485 | 512 | 542 | 145 | 137 | 130 |
| 11 | 595 | 631 | 670 | 130 | 121 | 115 |
| 12 | 1108 | 1193 | 1282 | 367 | 350 | 337 |
| 13 | 660 | 702 | 746 | 206 | 196 | 187 |
| 14 | 703 | 748 | 796 | 158 | 148 | 140 |
| 15 | 816 | 871 | 930 | 202 | 190 | 180 |
| 16 | 912 | 976 | 1045 | 195 | 182 | 172 |
| 17 | 893 | 956 | 1023 | 176 | 164 | 155 |
| 18 | 1020 | 1097 | 1178 | 169 | 156 | 147 |

Estimated accuracy = $\pm 3\%$

APPENDIX III

GEL T1 AND T2 VALUES

Imaging field 0.28 Tesla

| Temperature (K) | | | | | | |
|-----------------|---------|------|------|---------|-----|-----|
| | 292 | 296 | 300 | 292 | 296 | 300 |
| Tube No | T1 (ms) | | | T2 (ms) | | |
| 1 | 158 | 165 | 172 | 50 | 47 | 45 |
| 2 | 237 | 249 | 261 | 70 | 66 | 63 |
| 3 | 238 | 249 | 261 | 106 | 102 | 100 |
| 4 | 357 | 377 | 396 | 52 | 48 | 45 |
| 5 | 359 | 378 | 397 | 91 | 86 | 81 |
| 6 | 361 | 379 | 398 | 146 | 140 | 136 |
| 7 | 482 | 510 | 537 | 93 | 86 | 81 |
| 8 | 485 | 512 | 539 | 131 | 124 | 118 |
| 9 | 609 | 646 | 684 | 114 | 106 | 100 |
| 10 | 611 | 648 | 685 | 153 | 143 | 136 |
| 11 | 739 | 787 | 836 | 135 | 125 | 118 |
| 12 | 1284 | 1388 | 1496 | 382 | 363 | 348 |
| 13 | 813 | 867 | 921 | 217 | 205 | 195 |
| 14 | 860 | 919 | 979 | 164 | 153 | 144 |
| 15 | 983 | 1054 | 1126 | 210 | 196 | 186 |
| 16 | 1085 | 1167 | 1251 | 200 | 187 | 176 |
| 17 | 1065 | 1145 | 1227 | 181 | 168 | 158 |
| 18 | 1197 | 1293 | 1391 | 172 | 159 | 149 |

Estimated accuracy = $\pm 3\%$

APPENDIX III

GEL T1 AND T2 VALUES

Imaging field 0.50 Tesla

| Temperature (K) | | | | | | |
|-----------------|---------|------|------|---------|-----|-----|
| | 292 | 296 | 300 | 292 | 296 | 300 |
| Tube No | T1 (ms) | | | T2 (ms) | | |
| 1 | 188 | 202 | 217 | 52 | 49 | 47 |
| 2 | 282 | 304 | 326 | 72 | 69 | 66 |
| 3 | 281 | 304 | 326 | 112 | 109 | 107 |
| 4 | 422 | 456 | 490 | 53 | 49 | 46 |
| 5 | 422 | 456 | 490 | 94 | 88 | 84 |
| 6 | 422 | 456 | 491 | 153 | 148 | 144 |
| 7 | 562 | 608 | 656 | 95 | 88 | 83 |
| 8 | 562 | 609 | 657 | 136 | 128 | 122 |
| 9 | 702 | 762 | 823 | 116 | 108 | 102 |
| 10 | 703 | 763 | 824 | 157 | 147 | 140 |
| 11 | 843 | 917 | 993 | 137 | 127 | 120 |
| 12 | 1395 | 1531 | 1674 | 389 | 370 | 356 |
| 13 | 918 | 1000 | 1084 | 222 | 210 | 201 |
| 14 | 970 | 1058 | 1148 | 166 | 155 | 147 |
| 15 | 1097 | 1198 | 1303 | 213 | 200 | 189 |
| 16 | 1202 | 1315 | 1432 | 203 | 189 | 178 |
| 17 | 1183 | 1294 | 1409 | 183 | 170 | 160 |
| 18 | 1318 | 1444 | 1576 | 173 | 160 | 151 |

Estimated accuracy = $\pm 3\%$

APPENDIX III
GEL T1 AND T2 VALUES

Imaging field 1.00 Tesla

| Temperature (K) | | | | | | |
|-----------------|---------|------|------|---------|-----|-----|
| | 292 | 296 | 300 | 292 | 296 | 300 |
| Tube No | T1 (ms) | | | T2 (ms) | | |
| 1 | 198 | 219 | 241 | 52 | 50 | 48 |
| 2 | 296 | 327 | 360 | 73 | 69 | 67 |
| 3 | 294 | 326 | 359 | 113 | 111 | 109 |
| 4 | 449 | 496 | 546 | 53 | 49 | 46 |
| 5 | 443 | 490 | 540 | 94 | 89 | 85 |
| 6 | 441 | 488 | 538 | 154 | 150 | 147 |
| 7 | 592 | 654 | 721 | 95 | 89 | 84 |
| 8 | 589 | 651 | 718 | 136 | 129 | 124 |
| 9 | 738 | 816 | 900 | 116 | 108 | 102 |
| 10 | 734 | 813 | 896 | 157 | 149 | 142 |
| 11 | 883 | 978 | 1079 | 137 | 128 | 121 |
| 12 | 1433 | 1591 | 1760 | 390 | 373 | 358 |
| 13 | 954 | 1057 | 1167 | 223 | 212 | 203 |
| 14 | 1013 | 1122 | 1238 | 167 | 156 | 147 |
| 15 | 1139 | 1263 | 1395 | 214 | 201 | 190 |
| 16 | 1249 | 1385 | 1530 | 203 | 190 | 179 |
| 17 | 1233 | 1366 | 1508 | 183 | 171 | 161 |
| 18 | 1375 | 1524 | 1683 | 174 | 161 | 151 |

Estimated accuracy = $\pm 3\%$

APPENDIX III
GEL T1 AND T2 VALUES

Imaging field 1.50 Tesla

| Temperature (K) | | | | | | |
|-----------------|---------|------|------|---------|-----|-----|
| | 292 | 296 | 300 | 292 | 296 | 300 |
| Tube No | T1 (ms) | | | T2 (ms) | | |
| 1 | 199 | 221 | 245 | 52 | 50 | 48 |
| 2 | 298 | 331 | 367 | 73 | 70 | 67 |
| 3 | 295 | 329 | 365 | 113 | 111 | 110 |
| 4 | 455 | 506 | 560 | 53 | 49 | 46 |
| 5 | 446 | 496 | 550 | 94 | 89 | 85 |
| 6 | 442 | 492 | 546 | 154 | 151 | 147 |
| 7 | 597 | 664 | 736 | 95 | 89 | 84 |
| 8 | 592 | 659 | 731 | 136 | 129 | 124 |
| 9 | 745 | 828 | 918 | 116 | 108 | 102 |
| 10 | 739 | 822 | 911 | 157 | 149 | 142 |
| 11 | 892 | 992 | 1099 | 137 | 128 | 121 |
| 12 | 1440 | 1603 | 1779 | 390 | 373 | 359 |
| 13 | 959 | 1068 | 1184 | 223 | 212 | 203 |
| 14 | 1023 | 1137 | 1260 | 167 | 156 | 148 |
| 15 | 1149 | 1278 | 1416 | 214 | 201 | 191 |
| 16 | 1261 | 1403 | 1555 | 204 | 190 | 180 |
| 17 | 1246 | 1385 | 1535 | 183 | 171 | 161 |
| 18 | 1392 | 1547 | 1714 | 174 | 161 | 151 |

Estimated accuracy = $\pm 3\%$

APPENDIX III

GEL T1 AND T2 VALUES

Imaging field 2.00 Tesla

| Temperature (K) | | | | | | |
|-----------------|---------|------|------|---------|-----|-----|
| | 292 | 296 | 300 | 292 | 296 | 300 |
| Tube No | T1 (ms) | | | T2 (ms) | | |
| 1 | 199 | 222 | 247 | 52 | 50 | 48 |
| 2 | 298 | 333 | 370 | 73 | 70 | 67 |
| 3 | 296 | 330 | 367 | 113 | 111 | 110 |
| 4 | 459 | 511 | 566 | 53 | 49 | 46 |
| 5 | 448 | 499 | 554 | 94 | 89 | 85 |
| 6 | 443 | 494 | 549 | 154 | 151 | 148 |
| 7 | 600 | 669 | 743 | 95 | 89 | 84 |
| 8 | 594 | 662 | 736 | 136 | 129 | 124 |
| 9 | 749 | 834 | 926 | 116 | 109 | 103 |
| 10 | 741 | 826 | 918 | 157 | 149 | 142 |
| 11 | 897 | 999 | 1109 | 137 | 128 | 121 |
| 12 | 1444 | 1609 | 1787 | 390 | 373 | 359 |
| 13 | 962 | 1073 | 1191 | 223 | 212 | 203 |
| 14 | 1028 | 1145 | 1271 | 167 | 156 | 148 |
| 15 | 1154 | 1285 | 1427 | 214 | 201 | 191 |
| 16 | 1268 | 1412 | 1567 | 204 | 190 | 180 |
| 17 | 1253 | 1395 | 1548 | 183 | 171 | 161 |
| 18 | 1403 | 1561 | 1731 | 174 | 161 | 151 |

Estimated accuracy = $\pm 3\%$

APPENDIX III
GEL T1 AND T2 VALUES

Imaging field 3.00 Tesla

| Temperature (K) | | | | | | |
|-----------------|---------|------|------|---------|-----|-----|
| | 292 | 296 | 300 | 292 | 296 | 300 |
| Tube No | T1 (ms) | | | T2 (ms) | | |
| 1 | 200 | 223 | 249 | 52 | 50 | 48 |
| 2 | 299 | 334 | 372 | 73 | 70 | 67 |
| 3 | 296 | 331 | 368 | 113 | 111 | 110 |
| 4 | 463 | 516 | 574 | 53 | 49 | 46 |
| 5 | 450 | 502 | 559 | 94 | 89 | 85 |
| 6 | 444 | 496 | 552 | 154 | 151 | 148 |
| 7 | 604 | 674 | 750 | 95 | 89 | 84 |
| 8 | 596 | 666 | 741 | 136 | 129 | 124 |
| 9 | 754 | 841 | 935 | 116 | 109 | 103 |
| 10 | 745 | 831 | 924 | 157 | 149 | 142 |
| 11 | 903 | 1007 | 1120 | 137 | 128 | 121 |
| 12 | 1448 | 1615 | 1796 | 390 | 373 | 359 |
| 13 | 966 | 1078 | 1199 | 224 | 212 | 203 |
| 14 | 1034 | 1153 | 1282 | 167 | 156 | 148 |
| 15 | 1160 | 1293 | 1437 | 214 | 201 | 191 |
| 16 | 1276 | 1422 | 1581 | 204 | 190 | 180 |
| 17 | 1262 | 1407 | 1563 | 184 | 171 | 161 |
| 18 | 1415 | 1576 | 1750 | 174 | 161 | 151 |

Estimated accuracy = $\pm 3\%$